



Disposable Surgical Face Mask







Technical Specification

NAME

Disposable Surgical Face Mask

SPECICATION

Flat Ear loop 17.5 cm*9.5 cm

CERTIFICATION

CE & FDA

STANDARD EN 14683:2019 TYPE IIR & YY 0469-2011

INTRUCTIONS

Place The Mask On The Lower Half Of The Face, And The Ear Hanging Rope Outward. Press The Bridge Of The Nose To Ensure Full Fit. Pull The Lower End Of The Mask Down To Cover The Lower Jaw. Hook The Ear Rope To The Back Of The Ear. The Mask Completely Covers The Lower Half Of The Face And Obtains The Maximum Area Of Protection.

APLICATION SCOPE

It Is Suitable For The Health Care Of The Wearer In The General Medical Environment.

STRUCTURE AND COMPOSITION

This Product Is Composed Of Three Materials: Non-Woven Fabric, Nose Strip And Elastic Band. The Ear Strap Is Made Of Elastic Band, Which Is Made Of Non-Woven Fabric With An Elastic Band Inside. The Material Of The Nose Strip Is A Metal Strip, Which Is Covered With Fine Galvanized Iron Wire Material.

MATERIAL

Three Layers 3 ply

Inner Layer Non-Woven Fabric

Middle Layer

Ultra-Fine Polypropylene

Fiber Melt Player Febrie

Fiber Melt- Blown Fabric

Outer Layer

Non-Woven Fabric Or Ultra-Thin
Polypropylene Melt- Blown Fabric



